

Principles of targeted sequencing analysis (Mutation and Copy number alternations)

Terms

Mutation

permanent inheritable change in DNA base sequence

Variants/Mutation/Polymorphism

Variant: broadest term

Polymorphism : common ($>1\%$) variants in populations

Mutation: rare ($<1\%$) variants in populations

Driver mutation

Genetic change that give selective advantage to clones during cancer development

Passenger mutations

Genetic change that do not confer any selective advantage in cancer development

Germline mutations

Mutations in germ cells that are inherited

Somatic mutations

Mutations in non-germ cells that are acquired as opposed to inherited

Terms

Base substitutions / insertion and deletions / structural variations

Base substitutions: A type of mutation in which one base is replaced by another

Insertions and deletions (Indels): A type of mutation that arises from the insertion or deletion of one or more nucleotides

Structural variations: Large-scale genomic changes (typically >1kb) such as deletions, tandem duplications, amplifications, inversions and translocations.

Copy number variation/alternations

: Phenomenon in which sections of the genome (>1kb) are repeated and the number of repeats in the genome varies between individuals

: Difference between variants/alternations similar to polymorphism/mutations

→ further classified as amplification (> 3 copies), gain, neutral (2), loss, deletion (0)


Loss of heterozygosity (LOH)

Loss of heterozygosity is defined as the loss of one parent's contribution to the cell and can be caused by deletion, or others (not equal to loss)

Allelic imbalance (AI)

Phenomenon where the two alleles of a given gene are differentially expressed

Classification of mutations



Principle of mutation

Rare event == Likely to be more pathogenic

Polymorphism: different between populations (ExAC, gnomAD east asian)

(Procto)oncogene: with hotspot (ex BRAF V600E residue)

Tumor suppressor gene: without hotspot (ex TP53)

Two-hit hypothesis (The Knudson hypothesis): two hits required for cancer in TSG

Possible normal cell contamination in tumor (more likely)

Possible tumor cell contamination in normal (less likely)

Normal sample is needed for genotype, which could be revealed in low depth

VAFs (Variants allele fraction): read counts supporting variants / total read counts

VAFs can be affected by tumor purity, clonal evolution, CNV status


CCF (cancer cell fraction) inferred by VAFs after adjusting purity and CNV

Clonal (all/most of tumor cells have mutation), Subclonal (Not all)

FFPE (False C>T at low VAFs), cfDNA (Very low VAFs)

Principle of mutation


Recurrent gene mutations in study or in DB== Likely driver mutation
But by not only simple mutation counts (gene length, methylation, ...)



MutSig builds a model of the background mutation processes (BMR) that were at work during formation of the tumors, and it analyzes the mutations of each gene to identify genes that were mutated more often than expected by chance, given the background model.

MutSigCV (CV for 'covariate') improves the BMR estimation by pooling data from 'neighbor' genes with similar genomic properties such as DNA replication time, chromatin state (open/closed), and general level of transcription activity.

Principle of mutation



Tumor purity: fraction of cancer cells within sample

Read depth: read count at specific locus

VAFs: read counts supporting variants / total read counts

If tumor purity 100% and no CNAs, VAF will be 50%




Clinical implications

Table 1 | Genomic assays and the clinical decisions they support

Omic feature	Time of aberration discovery	Clinical setting	Analyte	Analysis	Treatment implication	Time of implementation of first targeted therapy
<i>BCR-ABL</i> translocation	1973 ⁴¹	Chronic myeloid leukaemia, acute lymphoblastic leukaemia	DNA, RNA	PCR (translocation)	Imatinib, dasatinib, nilotinib	2001 ⁴²
<i>BCR-ABL</i> imatinib resistance mutation	2002 ¹¹⁹	Chronic myeloid leukaemia	DNA	Sequencing (point mutation)	Dasatinib, nilotinib	2006 ¹²⁰
<i>PML-RAR</i> translocation	1990 ¹²¹	Acute promyelocytic leukaemia	DNA	PCR (translocation)	All <i>trans</i> retinoic acid	1987 ¹²²
<i>KIT</i> mutation	2001 ¹²³	Gastrointestinal stromal tumour	DNA	Sequencing	Imatinib	2001 ¹²³
<i>HER2</i> amplification	1985 ¹²⁴	Breast cancer, gastric cancer	Tissue	Fluorescence <i>in situ</i> hybridization, immunohistochemistry	Trastuzumab, lapatinib, pertuzumab	2001 ¹²⁵
Oestrogen or progesterone receptors	1896 ¹²⁶	Breast cancer	Tissue	Immunohistochemistry	Multiple hormonal-based therapies (for example, aromatase inhibitors)	1896 (oophorectomy) ¹²⁶
<i>EGFR</i> mutation	2004 ⁷⁶	Non-small-cell lung cancer	DNA	Sequencing or PCR (point mutations)	Erlotinib, gefitinib	2003 ¹²⁷
<i>ALK</i> fusion gene	2007 ⁴⁴	Non-small-cell lung cancer	DNA	Fluorescence <i>in situ</i> hybridization (break apart)	Crizotinib	2010 ⁴³
<i>BRAF</i> mutation	2001 ¹²⁸	Melanoma	DNA	Sequencing or PCR (point mutations)	Vemurafenib	2011 ¹²⁹
<i>KRAS</i> mutation	1987 ¹³⁰	Colorectal cancer	DNA	Sequencing or PCR (point mutations)	Withhold cetuximab or panitumumab	2006 ¹³¹

Mutations / CNAs / Fusions : all events might be clinically applicable


Cancer genetics (PCAWWG paper)




Most of cases, but not all cases have mutations in DNA

Mutational landscape differ by disease types as well as tumors of the same type

Cancer genetics (PCAWWG paper)



Cancer genetics (PCAWG paper)



Two-hit hypothesis (The Knudson hypothesis): two hits required for cancer in TSG
 TSG genes altered as biallelic in most cases

Article

Pan-cancer analysis of whole genomes

<https://doi.org/10.1038/s41586-020-1969-6> The ICGC/TCGA Pan-Cancer Analysis of Whole Genomes Consortium

Further reading (Nature, 2020)

수천 cases, pan-cancer, WGS, ICGC/PCAWG

Sequencing of DNA




Figure 1 | Sequencing-based technologies in biomedicine.

Targeted/WES/WGS

	Targeted	WES	WGS
Application	Mutation, CNAs (MSI, TMB)	Mutation, CNAs, Signature, MSI, TMB Purity, Telomere	WES applications + SVs (chromothripsis, translocation, fusion)
Mutation count	수~수십	수십~수천	수천~수십만
Target region	Various Less than 3Mb	3-6Mb	1Gb~3Gb
Sequencing	Various Less than 3Gb	About 10Gb	About 30Gb
Indication	cfDNA, small FFPE, Validation	Research standard Frozen/FFPE	Future standard Frozen Difficult to analysis and interpret
Kit	Ion AmpliSeq Celemics Kit SureSelectXT Target	Sureselect Exon V6	TruSeq DNA Nano
Cost	30만원후	80만원후	250만원후

Panel sequencing

Panel?

Which gene, mutation/CNAs, Hotspots/coding sequence

Pros

1. Focused on Clinically applicable information
2. Accuracy due to high read depth
3. High sensitivity detecting low VAF mutations
4. Small amount DNA available (10ng)
5. FFPE/Liquid biopsy available
6. Low cost
7. Faster Sequencing and Analysis


Panel sequencing only available (small sample, small FFPE, cfDNA,...)

Cons

mutation/CNAs only, outside panel unknown

Panel sequencing

A




Ampli-seq
(Thermo Fischer)

Hyb-seq
(Agilent, Celemics)


B

Amplicon-based assay
(Tiled amplicon approach)



C


Hybridization capture-based assay



Panel sequencing

Sample source	PCR based = Hybrid capture based (Tissue, FFPE, liquid biopsy)	
Input DNA	PCR based	> Hybrid capture based
Panel size	PCR based	< Hybrid capture based
Data generation time	PCR based	> Hybrid capture based
Homopolymer analysis	PCR based	< Hybrid capture based
Application	PCR based (Mutation, CNA)	Hybrid capture based (Mutation, CNA, MSI test)
	More prone to amplification error (?)	High off-target (Large CNAs) (but more cose)

Low sample input amount = More amplification, high duplication rate = low read depth



* 비유전성 고형암 필수 유전자

Mutation annotation and prioritization


Step1: Classify mutations into germline and somatic mutations

Step2: Annotate **somatic mutations** (add information)
using webannovar <http://wannovar.wglab.org/>

Step3: Filter less likely somatic mutations (Filter subject to research objective)
Population minor allele frequency > 0.1% (ExAC, gnomAD)
Read counts < 20 (Normal or tumor)
VAFs < 2-5%

Step4: Prioritize somatic mutations (to find putative driver mutations)
COSMIC, cbiportal, oncoKB, Polyphen-2, ClinVar,...

cfDNA/ctDNA



cfDNA


=cell-free DNA

: Degraded DNA
fragments released to
the blood plasma
: Including ctDNA



ctDNA

=circulating tumor DNA
:Tumor-derived
fragmented DNA in the
bloodstream

Ampliseq (1): Library preparation



Ampliseq (2): Templating and sequencing





Mean 150-200bp

Total: ~1Gb

Loading with Ion Chef™

Manual loading



Key point for successful DNA sequencing study

DNA amount, Quality (DIN), Purity, Sample condition (FFPE, Thawing)

Park WY, Sci Rep

Agilent Tapestation DIN for FFPE







Table 1. The DIN and the sequencing quality criteria obtained for the six samples shown in Figure 3.

Sample	DIN	On-target rate % > 70 %	10x Coverage rate % > 90 %	Deduplication rate %
1	4.1	79.7	98.6	66.0
2	5.2	78.0	99.0	74.2
3	5.0	80.7	98.1	87.1
4	2.2	23.8	94.0	65.0
5	2.9	51.5	97.7	95.7
6	1.9	47.4	96.6	97.1



Same amount, different depth

False negative (esp. in low VAF) (DIN, conta) or false positive (Amp error)